


We have presence in all over India.

The markers are a symbolic representation only.



Meditech Electronics

CONTACT INFO

J-81, Sector-4, DSIIIDC,
Bawana Indl. Area, Delhi-110039
Ph : 011-27762484, +91-9868165949,
For Sales : 9811165949, 9210333574,
9811225949
E-mail : sbnain@yahoo.com;

JOIN US AT:

sbnain@gmail.com
 sbnain@yahoo.com
 sbnain@hotmail.com
 sbnain2512
 http://meditech-electronics.blogspot.in/

WE ALSO AVAILABLE AT :

www.facebook.com/ledotlight
 9868165949
info.9868165949@gmail.com
 QQ: 1517286642

Note : Due to upgrade the product, we right to make changes on product design, Colour & specification without prior notice
Authorised Dealers / Distributors



COMMITTED TO QUALITY
COMMITTED TO YOU



www.meditechelectronics.com | www.ledotlight.com

MEDITECH ELECTRONICS LED TECHNOLOGY



Model : ME601

Full Spectrum Colour Rendition : Reduce Eye Strain by making the light do the work Instead of the eyes

Efficiency and long life : The high efficiency achieves more light with less energy.

Environment Friendly : LED uses less energy to produced large beam. LEDs are mercury free, do not contain harmful chemical.

Advanced Technology to achieve depth of illumination, colour temperature and colour indication index of harmony and unity. Imported and High Performance LED and Optics for Energy Saving and Superior Quality View


Special LED to give out Direct Illumination for Deep Cavity Operation



Focusable






Shadow Less



Controlled Colour Temperature



Memory Function



UNIQUE DESIGN, NATURAL & COMFORTABLE LIGHT



Model : ME401H


Technical Data

Intensity	1,35,000 Lux
Spot Diameter	100 - 150 mm
Depth Illumination	10 - 15 cm
Colour Temp.	3800-5200 K
Colour Rendering Index	93 Ra
Focusing	Manual
LED Average Life	>50000 Hrs.
Number of LED	84
Diameter of Light	690mm
Intensity Control	10~100 Digital
Power Supply	220V/50Hz AC
Power Consumption	55 Watt
Special Function	Memory Function

Optional Remote



NO HEAT
COOL LIGHT




Model : ME301H

Technical Data

Intensity	1,10,000 Lux
Spot Diameter	100 - 150 mm
Depth Illumination	10 - 15 cm
Colour Temp.	3800-5200 K
Colour Rendering Index	94 Ra
Focusing	Manual
LED Average Life	>50000 Hrs.
Number of LED	63
Diameter of Light	690mm
Intensity Control	10~100 Digital
Power Supply	220V/50Hz AC
Power Consumption	50 Watt
Special Function	Memory Function

Optional Remote



EXCELLENT ILLUMINANCE,
QUALITY ASSURANCE



Model : ME621


Technical Data

Intensity	1,40,000 Lux
Spot Diameter	100 - 150 mm
Depth Illumination	10 - 15 cm
Colour Temp.	3800-5200 K
Colour Rendering Index	91 Ra
Focusing	Manual
LED Average Life	>50000 Hrs.
Number of LED	98
Diameter of Light	680mm
Intensity Control	10~100 Digital
Power Supply	220V/50Hz AC
Power Consumption	65 Watt
Special Function	Memory

Optional : Remote Wireless Wall Control



EXCELLENT IN DESIGN
AND PERFORMANCE



Model : ME306

Technical Data	
Intensity	1,50,000 Lux
Spot Diameter	100 - 150 mm
Depth Illumination	10 - 15 cm
Colour Temp.	3800-5200 K
Colour Rendering Index	91 Ra
Focusing	Manual
LED Average Life	>50000 Hrs.
Number of LED	111
Diameter of Light	620mm
Intensity Control	10~100 Digital
Power Supply	220V/50Hz AC
Power Consumption	85 Watt
Special Function	Memory
Optional Remote	




SLEEK DESIGN,
HIGH PERFORMANCE




Model : ME401

Technical Data	
Intensity	1,25,000 Lux
Spot Diameter	100 - 150 mm
Depth Illumination	10 - 15 cm
Colour Temp.	3800-5200 K
Colour Rendering Index	91 Ra
Focusing	Manual
LED Average Life	>50000 Hrs.
Number of LED	76
Diameter of Light	600mm
Intensity Control	10~100 Digital
Power Supply	220V/50Hz AC
Power Consumption	50 Watt
Special Function	Memory
Optional Remote	



SIMPLE, CLEAR AND
AFFORDABLE LIGHT



Model : ME301

Technical Data

Intensity	1,00,000 Lux
Spot Diameter	100 - 150 mm
Depth Illumination	10 - 15 cm
Colour Temp.	3800-5200 K
Colour Rendering Index	91 Ra
Focusing	Manual
LED Average Life	>50000 Hrs.
Number of LED	57
Diameter of Light	600 mm
Intensity Control	10~100 Digital
Power Supply	220V/50Hz AC
Power Consumption	45 Watt
Special Function	Memory

Optional Remote



Meditech Electronics

07

SLEEK DESIGN,
HIGH PERFORMANCE



Model : ME403

Technical Data

Intensity	1,25,000 Lux
Intensity Satellite	1,00,000 Lux
Spot Diameter	100 - 150 mm
Depth Illumination	10 - 15 cm
Colour Temp.	3800-5200 K
Colour Rendering Index	91 Ra
Focusing	Manual
LED Average Life	>50000 Hrs.
Number of LED	76 + 57
Diameter of Light	580mm X 2
Intensity Control	10~100 Digital
Power Supply	220V/50Hz AC
Power Consumption	40 + 50 Watt
Special Function	Memory

Optional Remote



Meditech Electronics

08

INNOVATIVE DESIGN , MEET ALL
THE NEED OF MODERN SURGERY.



Model : ME1000E

Electric Operation Table

Technical Specification

Length of Table	1910 mm
Width of Table Top	522 mm
Height (Min-Max.)	750-1000 mm
Trendelenburg (Rev.)	25°
Lateral Tilt	±20°
Back Section	±60°
Head Section	±60°
Leg Section	0°-90°

Suitable for :

- General Surgery
- Vascular Surgery
- Cardiac Surgery
- Neurology Surgery
- Urology Surgery
- Gynecology
- Laparoscopy
- Plastic Surgery

Standard Accessories

- Anesthetic Screen
- Shoulder Supports with pad
- Arm Boards with pad
- Lateral support with pad
- Knee Crutches
- Water Proof rubber Mattress



INNOVATIVE DESIGN , MEET ALL
THE NEED OF MODERN SURGERY.



Model : ME800H

C-Arm Compatible Hydraulic Operation Table

Technical Specification

Length of Table	1910 mm
Width of Table Top	522 mm
Height (Min-Max.)	750-1000 mm
Trendelenburg (Rev.)	25°
Lateral Tilt	±20°
Back Section	±60°
Head Section	±60°
Leg Section	0°-90°


Suitable for :

- General Surgery
- Vascular Surgery
- Cardiac Surgery
- Neurology Surgery
- Urology Surgery
- Gynecology
- Laparoscopy
- Plastic Surgery


Standard Accessories

- Anesthetic Screen
- Shoulder Supports with pad
- Arm Boards with pad
- Lateral support with pad
- Knee Crutches
- Water Proof rubber Mattress







Normal Position




General Surgery




Trendelenburg




Lateral Tilt Left




Lower Back Section UP




Upper Back Section Up




Head Section Up




Leg Section Up




General Surgery




Rev. Trendelenburg




Lateral Tilt Right




Lower Back Section Down



Upper Back Section Down



Head Section Down



Leg Section Down

Meditech Electronics

11



Radiolucent Operating Table with Stand



Radiolucent Operating Table with Stand



Femur Nailing Support



Neuro Sujita Head Fixator for Sitting / Prone Positions



Head Rest device for prone, Spine & Lateral Position



Spinal Bridge for prone position



Neuro Skull Fixator Spine Position



Steinmen Pin Holding Clamp



Finger GripvDevice



Tibia Nailing

STANDARD ACCESSORIES



Lithotomy Crutches



Shoulder Support



Crutches Holder Steel Clamp



Lateral Support



Arm Rest



Water Proof Antistatic Rubber Mattress

Meditech Electronics

12

INNOVATIVE DESIGN , MEET ALL
THE NEED OF MODERN SURGERY.



Model : ME700G-E

Gynecology Operating Table

Awesome Design

Durable, lightweight and easy to use hand control with ergonomic shape featuring easy-to read hand control pictograms.

Fully Laminated Mattress

MOULDED WATERPROOF SEALED BELLOW protects against infiltration of body fluids.


Easy To Use

Minimum Height upto 60 cm provides easy access for lady surgeon.

High Quality Mattress

A high quality 3 inches mattress enhances the adaptive capacity of the mattress.

In terms of quality, reliability and functionality Operating Table is being acknowledged as the ideal solution for Gynaecology and Obstetrics. Due to its modular design, can be flexibly adapted to your needs and configured to suit your exact surgical requirements.



Just a few to begin with. More in the pipeline! Rest assured assured, even without them, we never compromise on quality!